# Tempcontrol 37-2 digital

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Knowledge of this manual is required for the operation of the device. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning the safe operation of the device.

Design and specifications are subject to change without notice. The manual is not covered by an update service.

Unless expressly authorized, forwarding and duplication of this document, and the utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

All product names mentioned herein may be the trademarks or registered trademarks of their respective companies and "TM" and "®" are not mentioned in each case in this manual.

© 2002

This symbol is a warning which must be observed under all circumstances!

This symbol is a warning which indicates information for the safety of the user and for the avoidance of harm to the devices.

This symbol is a warning which indicates to disconnect the instrument from the line!


**Before connecting to the line, the selected line voltage at the device is to be checked and to be changed if necessary! Put in the correct fuse!**

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




rear view

## Description

- The temperature regulator Tempcontrol 37-2 digital is designed for heating and controlling of special microscope stages, plates etc. (see below). The heating comes from the lost heat of transistors.
- The two heating channels can be activated separately and have a control range from 3°C over room temperature up to approx. 60°C.
- The device has two separated nominal values. The real values are displayed on a digital multifunctional display (display in °C; resolution 0.1°C)
- The RS 232-interface enables the control of the device by a computer.

Heatable Universal Mounting Frames A-H, K-H and M-H

Heating Insert for stereo microscope Zeiss

Heating Workplates

Heating Frame

CTI-Controller

Heating Microscope Stage

Heating and Cooling Frame

Heating Insert P and M 06-96

Objective Heater

## Specifications


AC operating voltage	115 / 230 V, 50 ... 60 Hz
Rated power	200 VA max.
Fuse	115 V – T3.15 A/E, 230 V – T1.6 A/E
Protective low voltage	20 V DC (output)
Control range	3°C above RT up to 60°C
Control tolerance	0.1°C
CE-Norm	class of protection 1
Dimensions (WxHxD)	175 x 138 x 270 (in mm)
Weight (kg)	5.0

## Accessories (included)

PC software "Remote Console Program"	# 0420.317
PC interface cable + adapter	# 0420.314
Power cable	# 0420.001

### Heatable Universal Mounting Frames

M-H # 000000-1116-054


A-H # 000000-1116-055

K-H # 000000-1116-056

KH-L # 000000-1005-837


### Heatable Microscope Stage

# 000000-1116-048

### Heating Insert P

# 411861-9901-000


### Objective Heater

# 000000-1116-073


# 000000-1116-075

# 000000-1116-076


## Assembly

- Ideally the device is placed on the right next to the microscope. The back side must have at least 10 cm free space to ensure sufficient cooling. The cooling is supported by a temperature controlled fan, which is switched-on at a temperature of approx. 70°C (at the transformer).
- Plug-in the 8-pin connector of the unit to be controlled into the socket of channel 1 or 2 at the back and tighten it.

## Operation

- Switch-on the Tempcontrol 37-2 digital with the power switch at the front panel.
- The following description assumes that the display shows the **temperature real values**. After completion of the input the display jumps automatically back to the real value after approx. 8 seconds. The changes made are then saved.
- The device has two independent channels, which can be controlled by the buttons **1** ("CH1" left in the display) or **2** ("CH2" left in the display), respectively.
- **Select control algorithm parameters:**

Because of the different characteristics of the connectable units it is absolutely necessary to select the correct control parameter set for the PID loop control.

Therefore when a new unit has been connected the control parameter set of the regarding channel must be selected first:


Control parameter set:	Suitable for:
PA1 (standard choice)	All heatable components, except Heating Unit
PA2	Heating Unit

**3x push button **1** or **2**:** the display shows then the control parameter set that is used for this channel: **PA1** or **PA2**

Use buttons **+** or **-** to choose the correct control parameter set.

The display jumps back to the real value after 8 s.

- **Set nominal value:**  
 2x push button **1** or **2**: the display shows then the nominal value that is used for this channel: **10.0°C to 60.0°C**  
 Use buttons **+** or **-** to choose the correct control parameter set.  
 The display jumps back to the real value after 8 s.
- **Switch the heating on or off:**  
 1x push button **1** or **2**: the display shows then the status of the heating for this channel: **on or off**  
 Push button **+** to switch-on the heating (green LED at button **1** or **2** lights up)  
 Push button **-** to switch-off the heating (green LED at button **1** or **2** goes off)  
 The display jumps back to the real value after 8 s.
- **Display / Operation scheme:**

## Remote-Operation

Press the **E** button (red LED at the button **E** lights up) to switch the device into the remote control mode. In this mode the device accepts no commands from the buttons at the front panel. The device is now exclusively controlled by the PC program "Remote Console Program".

After transition into the remote operation mode the device continues to work with the last parameters set. To deactivate the remote control mode press **E** again (red LED at button **E** goes off). The device can now again be controlled by the buttons **1** or **2**.

## Connection to the PC and installation of the operating software

- The Tempcontrol 37-2 digital is connected to the PC using the RS-232 interface cable supplied.
- Windows™ 3.1 or higher is required to use the software.
- The installation of the software is described in the "readme" file (3.5" Diskette)

## Operation via PC

- All Tempcontrol 37-2 digital functions (except for the power switch) can be controlled using the PC program "Remote Console Program".
- Before the program can communicate with the controller the interface address must be adjusted. To do this, the correct interface address must be selected in the menu "Controller>Select one".
- Operation by PC is similar to direct operation using the front panel. The control buttons are activated using the mouse.

## Maintenance

- The device is service-free.
- The unit can be cleaned from the outside with a lightly moistened cloth, if necessary.
- No water may get into the unit!
- Do not use aggressive or scrubbing cleaners!